
Invacare® Roze™

en	Stand Assist Patient Lift User Manual	4
da	Stålift Brugsanvisning	48
de	Aufstehhilfe und Patientenlifter Gebrauchsanweisung.....	92
es	Grúa de bipedestación Manual del usuario.....	138
fr	Lève-personnes verticalisateur Manuel d'utilisation	182
nl	Stand Assist Patiëntlift Gebruiksaanwijzing	226
no	Personløfter med ståfunksjon Bruksanvisning	272

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

Yes, you can.®

© 2014 Invacare Corporation.

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. A spare parts catalog is available upon request.

Contents

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.

I General	5
1.1 Symbols	5
1.2 Intended Use	6
1.2.1 Contraindications	6
1.3 Included Items	7
1.4 Service Life	8
1.5 Warranty Information	8
2 Safety	9
2.1 General Guidelines	9
2.2 Operating Information	9
2.2.1 General	9
2.2.2 Positioning	10
2.2.3 Disposal	11
2.3 Radio Frequency Interference	11
2.4 Product Labeling	12
3 Setup	13
3.1 Safe Assembly	13
3.2 Assembling the Mast to the Base	14
3.3 Assembling the Knee Pad to the Mast Assembly	15
3.4 Assembling the Foot Plate to the Mast Assembly	15
3.5 Checking the Service Light	16
3.5.1 Resetting the Service Light	16
3.6 Attaching the Battery Charger Mounting Bracket to the Wall	17
3.7 Attaching the Battery Charger Mounting Bracket to the Mast Assembly	17

4 Usage	19
4.1 Introduction	19
4.2 Closing/Opening Legs	19
4.2.1 Closing/Opening Electric Legs	19
4.3 Locking/Unlocking the Rear Casters	20
4.4 Raising/Lowering the Lift	20
4.4.1 Raising/Lowering an Electric Lift	20
4.5 Activating a Mechanical Emergency Release	20
4.5.1 Primary Emergency Release	20
4.5.2 Secondary Emergency Release	21
4.6 Performing an Emergency Stop	21
4.7 Charging the Battery	21
4.7.1 Battery Indicator	22
4.7.2 Using a Power Cord to Charge the Battery	24
4.7.3 Using the Battery Charger to Charge the Battery	25
5 Lifting the Patient	26
5.1 Safe Lifting	26
5.2 Positioning the Stand Assist Lift	28
5.3 Lifting the Patient	29
5.4 Moving the Patient	31
6 Transferring the Patient	32
6.1 Safe Transfer	32
6.2 Transferring to a Commode	33
6.3 Transferring to a Wheelchair	34
6.4 Transferring to a Bed	34
7 Troubleshooting	36
7.1 Troubleshooting Table	36
8 Maintenance	38
8.1 Maintenance Safety Inspection Checklist	38
8.2 Service Interval	39
8.3 Cleaning the Sling and the Lift	39
8.4 Reuse	40

8.5	Detecting Wear and Damage	40
8.6	Lubricating the Lift.	40
8.7	Replacing the Electric Actuator.	41
8.8	Replacing the Leg Actuators	42
8.9	Adjusting the Knee Pad Height	43
9	Technical data.	44
9.1	Patient Lift.	44

I General

I.1 Symbols

Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING!

– Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION!

– Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

IMPORTANT

– Indicates a hazardous situation that could result in damage to property if it is not avoided.

Gives useful tips, recommendations and information for efficient, trouble-free use.

Date of manufacture.

This product complies with Directive 93/42/EEC concerning medical devices.

The launch date of this product is stated in the CE declaration of conformity.

Read Manual

Audible tone when battery low. Refer to 4.7 Charging the Battery, page 21.

Caster lock.

Open/close Legs

Raise/lower the boom

Emergency stop

Safe Working Load

Ensure clip is fully closed and sling loop is completely on the hook before lifting the patient.

Double Insulated, Class II equipment

Type B applied part

Recycle this product. Refer 2.2.3 Disposal, page 11.

I.2 Intended Use

WARNING!

Risk of Falling

The Invacare stand assist lift is **NOT** a transport device. It is intended to transfer an individual from one resting surface to another (such as a bed to a wheelchair). Invacare slings and patient lift accessories are specifically designed to be used in conjunction with Invacare patient lifts.

To avoid injury or falls:

- Patients who are transferred and positioned using the stand assist lift **MUST** be cooperative, coherent and have head and neck control. Otherwise, injury may occur.
- Patients who are transferred and positioned using the stand assist sling **MUST** be able to support a majority of their own weight. Otherwise, injury or damage may occur. Refer to the sling user manual for more information.

Stand assist lifts are battery-powered transfer devices, designed to be used in most of the common lifting situations, for example:

- Between the bed and a wheelchair
- To and from the toilet

The stand assist lift is only intended to lift patients up to the maximum weight limit stated in technical data.

Selecting the appropriate slings and accessories for each individual is important to assure safety when using a patient lift. See Invacare's sling and accessory user manuals for further information on those devices.

Invacare recommends that the patient be transferred to a shower chair or other means for bathing.

The stand assist lift can be turned (rotated) in place for transfers with limited floor space.

I.2.1 Contraindications

The stand assist patient lift is contraindicated for patients who:


- are not cooperative
- are not coherent
- do not have head and neck control

Some slings designed for the stand assist patient lift are contraindicated for patients who cannot support the majority of their weight. Refer to the sling user manual for more information.

1.3 Included Items

The items listed in the tables are included in your package. Slings are sold separately. The wall charger may be sold separately.

Item	Description	Quantity
A	Lifter	1
B	Battery	1
C	Charging Cord, EU	1
D	Charging Cord, UK	1
E	Hand Pendant (not shown)	1
F	Lifter User Manual (not shown)	1

The diagram shows a line drawing of a manual hydraulic lifter. Callout A points to the handle assembly at the top. Callout B points to the battery pack mounted on the side of the frame. Callout C points to a charging cord with a two-prong European-style plug. Callout D points to a charging cord with a three-prong British-style plug.

I.4 Service Life

WARNING!

Risk of injury

Maintenance **MUST** be performed only by qualified personnel.

Improper assembly may cause injury or damage.

- Regular maintenance of patient lifts and accessories is necessary to assure proper operation.
- **DO NOT** overtighten the mounting hardware. This will damage the mounting bracket.

The stand assist lift has an expected lifetime of 8 years when used in accordance with safety instructions, maintenance intervals and correct use, stated in this manual. The effective product lifetime can vary according to frequency and intensity of use.

Actuator Service Life

Number of lifts per day	Service life of actuator (in years)
1-2	10
3	9
4	6
5	5
6	4
7-9	3
10-13	2
14-27	1

I.5 Warranty Information

Terms and conditions of the warranty are part of the general terms and conditions particular to the individual countries in which this product is sold.

Contact information for your local Invacare office is located inside the back cover of this manual.

2 Safety

2.1 General Guidelines

WARNING!

– DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment - otherwise, injury or damage may occur.

ACCESSORIES WARNING

Invacare products are specifically designed and manufactured for use in conjunction with Invacare accessories.

Accessories designed by other manufacturers have not been tested by Invacare and are not recommended for use with Invacare products. In certain instances, the use of other manufacturer's slings may be possible.

– Contact your local Invacare office for more information about accessories.

NOTICE

– The information contained in this document is subject to change without notice.

Check all parts for shipping damage before using. In case of damage, DO NOT use the equipment. Contact the Dealer or Invacare representative for further instructions.

2.2 Operating Information

This section of the manual contains general safety information about your product. For specific safety information, refer to the appropriate section of the manual and procedures within that section. For instance, for safety information related to assembling the lift, refer to section 3 Setup, page 13.

2.2.1 General

WARNING!

Risk of Falling

DO NOT attempt any transfer without approval of the patient's physician, nurse or medical assistant. Thoroughly read the instructions in this user manual, observe a trained team of experts perform the lifting procedures and then perform the entire lift procedure several times with proper supervision and a capable individual acting as a patient.

- Use common sense in all lifts. Special care **MUST** BE taken with people with disabilities who cannot cooperate while being lifted.
- Use steering handle on the mast at all times to push or pull the patient lift.
- Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode).

WARNING!

Risk of Injury or Damage

The patient lift can be used indoors or outdoors. Some surfaces may cause the lift to be unstable and cause injury or damage.

The lift can be used in the shower or bath area, but excessive moisture will damage the lift and may cause injury.

- Avoid using the lift on an incline. Invacare recommends that the lift only be used on a flat surface.
- DO NOT roll casterbase over uneven surfaces that may cause the patient lift to tip over.
- DO NOT use the patient lift in the shower or bath or in any prolonged moisture environment.
- Ensure that the patient lift is wiped clean of any moisture after use.
- DO NOT store the lift in a damp area or in a damp condition. Refer to the storage temperature, humidity and pressure ranges listed in 9.1 Patient Lift, page 44.
- Periodically inspect all components of the patient lift for signs of corrosion. Replace all parts that are corroded or damaged.


2.2.2 Positioning

WARNING!

Risk of Injury

- ALWAYS be aware of the Lift Arms. Injury to the patient and/or assistant may occur.
- ALWAYS be aware of the Footrest, especially the patient’s position on the footrest. Injury to the patient and/or assistant may occur.

2.2.3 Disposal

**WARNING!****Environmental Hazard**

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EC.

Device contains lead acid batteries.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- DO NOT dispose of batteries in normal household waste. They MUST be taken to a proper disposal site. Contact your local waste management company for information.
- Please be environmentally responsible and recycle this product through your recycling facility at its end of life.

2.3 Radio Frequency Interference

**WARNING!**

- Most electronic equipment is influenced by Radio Frequency Interference (RFI). CAUTION should be exercised with regard to the use of portable communication equipment in the area around such equipment. If RFI causes erratic behavior, PUSH the Red Power Switch OFF IMMEDIATELY. DO NOT turn the Power Switch ON while transmission is in progress.

2.4 Product Labeling

⚠ WARNING Risk of injury or damage. Only use the Roze mast with this base.

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

3.1 Safe Assembly

**WARNING!****Risk of Injury**

Improper assembly may cause injury or damage.


- Assembly **MUST** be performed only by qualified personnel.
- Use only Invacare parts in the assembly of this patient lift. The lift components are manufactured to specifications that assure correct alignment of all parts for safe functional operation.
- **DO NOT** overtighten the mounting hardware. This will damage the mounting bracket.


There are no tools required to assemble the patient lift.

If there are any issues or questions during assembly, contact a local Invacare representative. Refer to the contact information in the back of this manual.


3.2 Assembling the Mast to the Base

DETAIL A

DETAIL B

Ⓐ Mast Assembly	Ⓕ Mast Assembly Cable
Ⓑ Base	Ⓖ Base Cable
Ⓒ Caster	Ⓗ Tube of the Mast Assembly
Ⓓ Top Knob	Ⓘ Hole in the Base
Ⓔ Bottom Knob	

The mast assembly may be removed from the base for storage or transporting. The mast assembly **MUST** be properly secured to the base assembly before use.

1. Put the base on a level surface so all casters make contact with the floor.
2. Lock the rear casters (Detail A).
3. Unscrew the top knob and remove it from the base.
4. Loosen the bottom knob, but leave it screwed into the base.
5. Connect the mast assembly cable to the base cable (Detail B).

WARNING!

Risk of injury

– DO NOT place hands between the mast and base during assembly.


CAUTION!

Risk of Damage.

– DO NOT crush the connected cables between the mast and base during assembly.

6. Put the tube of the mast assembly into the hole in the base.
7. While supporting the mast assembly, tighten the bottom plastic knob.
8. Screw the top plastic knob into the hole in the base.

3.3 Assembling the Knee Pad to the Mast Assembly


Ⓐ Knee Pad	Ⓒ Fix Block
Ⓑ Knee Pad Mounting Bracket	Ⓓ Quick Release Knob

The knee pad may be removed from the patient lift for storage or transporting. The knee pad **MUST** be properly secured to the mast assembly before use.

1. Push the knee pad down into position on the knee pad mounting bracket on the mast assembly.
2. Put the fix block in place on the knee pad mounting bracket.
3. Insert the quick release knob through the mounting bracket, knee pad and fix block.
4. Tighten the quick release knob.


3.4 Assembling the Foot Plate to the Mast Assembly

Ⓐ Foot Plate (shown in place)	Ⓒ Mast Assembly
Ⓑ Top Mounting Bolt	Ⓓ Bottom Mounting Bolt

1. Rotate the foot plate up and place the hooks on the foot plate onto the top mounting bolts in the mast assembly.
2. Rotate the foot plate down towards the mast assembly. The bottom of the foot plate will rest against the bottom bolts on the mast assembly.
3. Ensure that the foot plate is level and secure.

3.5 Checking the Service Light

 Each time the lift is assembled, and before using the lift, the service light should be checked.

1. Assemble the lift.
2. Examine the control box (A) or the hand pendant (B) to see if the service light (C) is flashing.
3. Refer to the table.


Service Light	Initial Assembly	Reassembly
Flashing	Reset the service light. Refer to 3.5.1 Resetting the Service Light, page 16.	The lift requires service. Contact your local Invacare dealer or representative for service.
Not Flashing	The lift is ready for use.	

3.5.1 Resetting the Service Light

This procedure should only be performed by a qualified technician.

This procedure should only be performed after performing 3.5 Checking the Service Light, page 16.

1. Locate the handset (A).
2. Press and hold the UP button (B) and DOWN button (C) at the same time for five seconds.


You will hear a sound when the service light has been reset.

The service light is now reset to indicate a service requirement after a 12 month period or 8000 cycles.

If the service light is required to be set for a period shorter than 12 months, a special handset and further instructions are required. Contact Invacare to order this handset and instructions.

3.6 Attaching the Battery Charger Mounting Bracket to the Wall


i Refer to your local regulations concerning proper mounting procedures.

1. Place the battery charger with mounting bracket **A** on the wall at the desired position.
2. With a pencil, mark the middle hole **B** position.
3. Measure down 6½ inches (16.5 cm) from the pencil mark and drill one mounting hole.
4. Install the bottom mounting screw **C** until there is an approximate 1/8-inch (3 mm) gap between the screw head and the wall.
5. Install the battery charger with mounting bracket onto the bottom mounting screw.
6. Drill the remaining two mounting holes.
7. Install the two remaining mounting screws **D** through the mounting bracket and into the wall. Tighten securely.
8. Plug the battery charger into the wall electrical outlet.

i ON LED should illuminate.

3.7 Attaching the Battery Charger Mounting Bracket to the Mast Assembly

i **A** = Battery and battery charger shown in place on mast assembly.

WARNING!

– If the battery charger is attached to the mast assembly and is plugged into a wall outlet, make sure to unplug the battery charger from the wall outlet **BEFORE** using the patient lift. Failure to unplug the battery charger may result in injury or damage.

1. Remove the two mylar pieces covering the mounting holes on the back of the mast assembly.
2. Insert the hook screws into mounting holes.
3. Tighten the hook screws.
4. Mount the battery charger mounting bracket on the hook screws.
5. Charge the battery. Refer to 4.7 Charging the Battery, page 21.

4 Usage

4.1 Introduction

The operation of the patient lift is an easy and safe procedure.

 Before using the lift with a patient, refer to the following procedures for safety information and instruction:

- 2.2 Operating Information, page 9
- 5.3 Lifting the Patient, page 29

4.2 Closing/Opening Legs

WARNING! **Risk of Injury**

The lift could tip and endanger the patient and assistants.


- The legs of the lift must be in the maximum open position for optimum stability and safety. If it is necessary to close the legs of the lift to maneuver the lift under a bed, close the legs of the lift only as long as it takes to position the lift over the patient and lift the patient off the surface of the bed. When the legs of the lift are no longer under the bed, return the legs of the lift to the maximum open position.

4.2.1 Closing/Opening Electric Legs

The pendant is used to open or close the legs of the base for stability when lifting a patient.

Refer to the safety information in 4.2 Closing/Opening Legs, page 19 before performing this procedure.


1. To close the legs, press and hold the legs closed button **A**.
2. To open the legs, press and hold the legs open button **B**.

The legs will stop moving when the button is released.

4.3 Locking/Unlocking the Rear Casters

- To Lock the Caster: Push tab down.
- To Unlock the Caster: Push tab up.

4.4 Raising/Lowering the Lift


WARNING!
Risk of Injury

The lift could tip and endanger the patient and assistants.
– Invacare does recommend that the rear casters be left unlocked during lifting procedures to allow the patient lift to stabilize itself when the patient is initially lifted from a chair, bed or any stationary object.

4.4.1 Raising/Lowering an Electric Lift

Refer to the safety information in 4.4 Raising/Lowering the Lift, page 20 before performing this procedure.

1. To raise the lift — Press and hold the UP **A** button to raise the boom and the patient.
2. To lower the lift — Press and hold the DOWN **B** button to lower the boom and the patient.


Release the button to stop raising or lowering the lift.

4.5 Activating a Mechanical Emergency Release


There are two types of mechanical emergency release — primary and secondary.

4.5.1 Primary Emergency Release

1. Insert a pen into the hole labeled Emergency Up **A** or Emergency Down **B** on the control box **C**.

4.5.2 Secondary Emergency Release

Ⓐ RED Emergency Grip

Ⓑ Lift Arms


It is recommended that the primary emergency release be used. The secondary emergency release is only a back-up to the primary emergency release.

In cases where the primary release is either not functioning or unreachable, a secondary emergency release may be used.

1. Pull up on the EMERGENCY grip Ⓐ and push down on the lift arms Ⓑ at the same time.

4.6 Performing an Emergency Stop

1. Press the RED button Ⓐ on the control box Ⓑ in to stop the boom and patient from raising or lowering.
2. To reset, rotate the emergency button clockwise.


4.7 Charging the Battery

Invacare recommends the battery be recharged daily to prolong battery life.

There are two different methods to charge the battery. One method uses a power cord that attaches to the control box, the other requires the battery to be mounted to the battery charger. Follow the appropriate procedure to charge the battery for your patient lift.

4.7.1 Battery Indicator

The battery indicator ① is located on the control box ② and on the hand pendant ③. The LEDs indicate the battery state:

Hand Pendant Battery Indicator ①	Control Box Battery Indicator ①	Battery State	Description
		Full Charge	The battery is OK — no need for charging (100–50%). The third LED is GREEN.
		Partial Charge	The battery needs to be charged (50–25%). The second LED is YELLOW.

Hand Pendant Battery Indicator [Ⓐ]	Control Box Battery Indicator [Ⓐ]	Battery State	Description
		Low Charge	The battery needs to be charged (Less than 25%). The horn beeps when a button is pressed. The first LED is YELLOW.
		Low Charge (LED blinking)	<p>The battery needs to be charged.</p> <p>Some of the functionality of the lift is lost and it is only possible to lower the boom.</p> <p>An audible alarm will sound (horn will beep) when battery is low. If the audible alarm sounds during a transfer, complete that transfer and then charge the battery.</p>

4.7.2 Using a Power Cord to Charge the Battery

A

A = Power cord connects here.


CAUTION!

– The emergency stop **MUST NOT** be activated - otherwise it will be impossible to charge the battery. While charging takes place, the patient lift cannot be used. **DO NOT** attempt to move the patient lift without unplugging the power cord from the wall outlet. **DO NOT** attempt to use the patient lift if the battery housing is damaged. Replace a damaged battery housing before any further use.

1. Attach the power cord to the control box.
2. Plug the power cord into a power outlet.
 -  The battery will charge in approximately 4 hours. Charging must be done in a room with good air ventilation.
3. Disconnect the power cord from the power outlet after the battery has been fully charged.

4.7.3 Using the Battery Charger to Charge the Battery

1. Lift up on the handle (A) on the back of the battery (B).
2. Lift the battery up and out away from the control box (C).

CAUTION!

Mounting the battery improperly may cause injury or damage.

- Make sure there is an audible click when mounting battery on the battery charger to confirm proper mounting.

3. Place the battery on the battery charger (D) as shown. Make sure there is an audible click.

The charge LED will illuminate. When charging is complete, charge LED will stop illuminating.

A battery needing to be fully recharged will take approximately four hours.

4. Lift up on the handle on the back of the battery.
5. Lift the battery up and out away from the battery charger.

CAUTION!

Mounting the battery improperly may cause injury or damage.

- Make sure there is an audible click when mounting battery on the control box to confirm proper mounting.

6. Reinstall the battery onto the control box as shown. Make sure there is an audible click.

The battery mounts to the control box and battery charger as shown.

5 Lifting the Patient

5.1 Safe Lifting

WARNING!

- DO NOT exceed the maximum weight limitation (SWL) of 200 kg (31 stone). Maximum patient weight is 200 kg (31 stone).
- DO NOT attempt to any transfer without approval of a Healthcare Professional.
- ALWAYS keep hands and fingers clear of moving parts to avoid injury.
- During transfer, with the patient suspended in the sling, DO NOT roll the base of the lift over any uneven surfaces that would cause the lift to become unstable.
- Use the steering handle on the mast assembly at all times to push or pull the stand assist lift.

WARNING!

Operating the Lift and Lifting the Patient

- Ensure equipment used with the stand assist lift is of adequate strength to lift the load (e.g. sling). Ensure that should the equipment malfunction that the person being lifted is not exposed to danger.
- Before positioning the legs of the stand assist lift around the patient, make sure that the patient's feet are out of the way of the foot plate, otherwise injury may occur.
- Adjustments for safety and comfort should be made before moving the patient. Patient's arms should be outside of the sling straps.

WARNING!

Operating the Lift and Lifting the Patient

- Before lifting a patient from a stationary object (wheelchair, commode or bed), slightly raise the patient off the stationary object and check that all sling attachments are secure. If any attachment is not correct, lower the patient and correct the problem, then raise the patient and check again.
- During transfer, with the patient suspended in a sling attached to the lift, DO NOT roll caster base over uneven surfaces that would create an imbalance of the patient lift and could cause the patient lift to tip over. Use steering handle on the mast at ALL times to push or pull the patient lift.
- Invacare recommends locking the rear swivel casters ONLY when positioning or removing the sling (stand assist or transfer) from around the patient.
- Invacare does NOT recommend locking of the rear casters of the patient lift when lifting an individual. Doing so could cause the lift to tip and endanger the patient and assistants. Invacare DOES recommend that the rear casters be left unlocked during lifting procedures to allow the patient lift to stabilize itself when the patient is initially lifted from a chair, bed or any stationary object.

**WARNING!****Using the Sling**

- Use the sling that is recommended by a Healthcare Professional for the comfort and safety of the individual that is being lifted.
- Individuals that use the stand assist sling **MUST** be able to support the majority of their own weight, otherwise injury may occur.
- Bleached, torn, cut, frayed or broken slings are unsafe and their use could result in injury or damage. Discard any sling that meets these criteria **IMMEDIATELY**.
- **DO NOT** alter slings.
- Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode).
- If the patient is in a wheelchair, activate the wheelchair brakes to prevent the chair from moving forwards or backwards.
- Consideration should be taken that prevents exposure to danger to the person being lifted, and that the person can be freed without injury.

**WARNING!****Stand Assist Sling**

- **DO NOT** use the stand assist sling in combination with the patient lift as a transport device. It is intended to transfer an individual from one resting surface to another (such as a bed to a wheelchair).
- Before lifting the patient, make sure the bottom edge of the stand assist sling is positioned on the lower back of the patient and the patient's arms are outside the stand assist sling. The belt **MUST** be snug, but comfortable on the patient, otherwise the patient can slide out of the sling during transfer possibly causing injury.
- The belt **MUST** be snug, but comfortable on the patient, otherwise the patient can slide out of the sling during transfer, possibly causing injury.

Transfer Stand Assist Sling

- Before lifting the patient, make sure the bottom edge of the transfer sling is at the base of the spine and the patient's arms are outside the transfer sling.
- **DO NOT** raise the patient to a full standing position while using the transfer sling, otherwise injury may occur.

Refer to the patient sling brochure for more information.

5.2 Positioning the Stand Assist Lift

Refer to the Safety section in this manual and review the information in 5.1 Safe Lifting, page 26 before proceeding further and observe all warnings indicated.


Before positioning the legs of the patient lift under a bed, make sure that the area is clear of any obstructions.

WARNING! **Risk of Injury**

The lift could tip and endanger the patient and assistants.

- The legs of the lift must be in the maximum open position for optimum stability and safety. If it is necessary to close the legs of the lift to maneuver the lift under a bed, close the legs of the lift only as long as it takes to position the lift over the patient and lift the patient off the surface of the bed. When the legs of the lift are no longer under the bed, return the legs of the lift to the maximum open position.

(A) UP Button	(D) Legs Close Button
(B) DOWN Button	(E) Hand Control
(C) Legs Open Button	

1. Press the legs open button on the hand control to open the legs of the stand assist lift to the maximum open position.
2. Position the stand assist lift using the mast handle.
3. Press the down arrow button on the hand control to lower the lift arms for easy attachment to the sling.


5.3 Lifting the Patient


WARNING!

- DO NOT exceed the maximum weight limitation (SWL) of 200 kg (31 stone). Maximum patient weight is 200 kg (31 stone).
- Individuals that use the stand assist sling MUST be able to support the majority of their own weight, otherwise injury may occur.
- DO NOT lock the rear casters of the stand assist lift when lifting an individual. Locking the rear casters could cause the stand assist lift to tip and endanger the patient and assistants. Unlocking the rear casters during lifting procedures allows the lift to stabilize itself while the patient is initially lifted from a wheelchair, bed of any stationary object.
- DO NOT move the patient if the sling is not properly attached to the attachment points on the stand assist lift.
- Make sure that the sling is properly attached BEFORE lifting the patient. If any attachments are not properly in place, correct the problem. When the sling is elevated a few inches/cm off the stationary surface and before moving the patient, check again to make sure the sling and all attachments are securely in place. If any problem is found, lower the patient back onto the stationary surface and correct the problem - otherwise, injury or damage may occur.
- Adjustments for safety and comfort should be made before moving the patient.


Detail A — Stand Assist Lift


Detail B — Sling Attachment

Detail C — Brakes

Detail D — Lifting the Patient

(A) Hand Grips	(F) Rear Caster
(B) Knee Pad	(G) Leg
(C) Foot Plate	(H) Brake

ⓓ Hooks	① Loop of Sling
ⓔ Clips	

 The patient **MUST** be in a seated position first. Use the head section of the bed to get the patient upright then move legs over the side of the bed.

1. Instruct the patient to hold onto the hand grips on both sides of the stand assist lift (Detail A).
2. Instruct the patient to lean back into the standing or transfer sling.

WARNING!

- Stand Assist Slings - Before lifting the patient, make sure the bottom edge of the sling is positioned on the lower back of the patient and the patient’s arms are outside the sling.
- Transfer Stand Assist Slings - Before lifting the patient, make sure the bottom edge of the sling is at the base of the spine and the patient’s arms are outside the sling.

3. Ensure the following:

- a. The patient’s knees are secure against the knee pad.
- b. The patient’s feet are properly positioned on the foot plate.
- c. Slings:
 - the bottom edge of the stand assist sling is positioned on the patient’s lower back.
 - the bottom edge of the transfer stand assist sling is at the base of the patient’s spine.
 - the patient’s arms are outside the sling.
 - the loops of the sling are completely on the hooks of the lift arms and the clip is fully closed (Detail B).
- d. The rear casters are unlocked.
- e. The legs are in the maximum open position.

WARNING!

- If transferring a patient from a wheelchair, the wheelchair brakes **MUST** be in the locked position before lowering the patient into the wheelchair. Otherwise, injury may occur.

4. If transferring from a wheelchair, engage the brakes on the wheelchair (Detail C).
5. Press the UP arrow button on the hand control to raise the patient above the surface (bed, wheelchair or commode). The patient’s weight will be fully supported by the stand assist lift (Detail D).

The lower center of gravity provides stability making the patient feel more secure and the lift easier to move.

5.4 Moving the Patient

WARNING!

- The legs of the stand assist lift **MUST** be in the maximum open position for optimum stability and safety. If the patient is in a sling and it becomes necessary to move through a narrow passage, close the legs of the stand assist lift only as long as it takes to move through the passage. When the stand assist lift is through the passage, return the legs to the maximum open position.
- **DO NOT**, during transfer of a patient suspended in the lift sling, roll caster base over uneven surfaces that would create an imbalance of the lift. This could cause the lift to tip over. Use the mast handle at all times to push or pull the lift.

1. Ensure the legs of the stand assist lift are in the maximum open position. If not, press the **OPEN LEGS** button on the hand control to move the legs to the maximum open position.
2. Move the stand assist lift away from the surface they were lifted from.
3. Slowly move the patient to the desired surface.

6 Transferring the Patient

6.1 Safe Transfer

WARNING!

- DO NOT exceed the maximum weight limitation (SWL) of 200 kg (31 stone). Maximum patient weight is 200 kg (31 stone).
- DO NOT attempt to transfer a patient without approval of a Healthcare Professional. Thoroughly read the instructions in this owner's manual, observe a trained team of experts performing the lifting procedures and then perform the entire lift procedure several times with proper supervision and a capable individual acting as a patient. Training can be provided. Contact your local Invacare office for more information.
- DO NOT move the patient if the sling is not properly connected to the attachment points of the stand assist lift. Check that the sling is properly connected to the attachment points prior to lifting a patient. If any attachments are not properly in place, correct the problem. When the sling is elevated a few inches/cm off the stationary surface and before moving the patient, check again to make sure that all sling attachments are secure. If any attachments are not properly in place, lower the patient back onto the stationary surface and correct this problem - otherwise, injury or damage may occur.
- Adjustments for safety and comfort should be made before moving the patient. The patient's arms should be outside the straps.

WARNING!

- Use the sling that is recommended by a Healthcare Professional for the comfort and safety of the individual that is being lifted.
- DO NOT lock the rear casters of the stand assist lift when lifting an individual. Locking the rear casters could cause the stand assist lift to tip and endanger the patient and assistants.
- The legs of the stand assist lift **MUST** be in the maximum open position for optimum stability and safety. If the patient is in a sling and it becomes necessary to move through a narrow passage, close the legs of the stand assist lift only as long as it takes to move through the passage. When the stand assist lift is through the passage, return the legs to the maximum open position. If it is necessary to close the legs to maneuver the stand assist lift under a bed, close the legs only as long as it takes to position the stand assist lift over the patient and lift the patient off the surface of the bed. When the legs of the stand assist lift are no longer under the bed, return the legs to the maximum open position.
- Be sure to check the sling attachments each time the sling is removed and replaced to ensure that it is properly attached before the patient is removed from a surface.

The use of one assistant is based on the evaluation of the health care professional for each individual case.

6.2 Transferring to a Commode

WARNING!


– Invacare recommends locking the rear swivel casters only when positioning or removing the sling from around the patient.

1. Lift the patient from the side of the bed.
2. Press the UP button on the hand control to elevate the patient high enough to clear the arms of the commode chair. Their weight will be supported by the stand assist lift.
3. Guide the patient onto the commode chair. This may require two assistants.
4. Press the down arrow button to lower the patient onto the commode chair.
5. Lock the rear swivel casters on the stand assist lift.


6. Depending on the style of sling, perform one of the following:
 - Stand Assist Sling - unhook the stand assist sling from the attachment points on the stand assist lift.
 - Transfer Stand Assist Sling -
 - a. Unhook the transfer stand assist sling from the bottom attachment points on the stand assist lift.
 - b. Lift up on the patient's legs and remove the thigh supports from underneath the patient.
 - c. If desired, unhook the transfer stand assist sling from the top attachment points on the stand assist lift.
-  The patient can remain in the upper portion of the transfer stand assist sling while using the commode.
7. Instruct or assist the patient in lifting their feet off the foot plate.
8. Remove the sling from around the patient.
9. Unlock the rear casters and pull the stand assist lift away from the commode.
10. When complete, recheck the sling for correct attachments.

6.3 Transferring to a Wheelchair

Detail A

Detail B

 **H** = Wheelchair Brake

1. Ensure the legs of the lift with the patient in the sling are in the open position. Press the legs open button until in maximum open position.
2. Move the wheelchair into position.

3. Engage the brakes of the wheelchair to prevent movement of the chair (Detail A).

WARNING!

– DO NOT place the patient in the wheelchair if the brakes are not engaged. The wheelchair brakes MUST be in a locked position before lowering the patient into the wheelchair for transport. Otherwise, injury may result.

4. Position the patient over the wheelchair (Detail “B”).
5. Press the down arrow button and lower the patient into the wheelchair.

WARNING!

– Lock the rear swivel casters ONLY when positioning or removing the sling (stand assist or transfer stand assist) from around the patient. Otherwise, injury may occur.

6. Lock the rear swivel casters.
7. Unhook the sling from all attachment points on the lift.
8. Instruct patient to lift their feet off the foot plate. Assist the patient if necessary.
9. Remove the sling from around the patient.
10. Unlock the rear casters and pull the stand assist lift away from the wheelchair.

6.4 Transferring to a Bed

The lower center of gravity provides stability making the patient feel more secure and the lift easier to move.

1. Position the patient as far over the bed as possible.

 If patient is being transferred from a surface that is lower than the bed, press the up arrow button to raise the patient above the surface of the bed. The patient should be elevated just high enough to clear the bed with their weight fully supported by the lift.

2. Press the down arrow button and lower the patient onto the bed.

**WARNING!**

– Invacare recommends locking the rear swivel casters **ONLY** when positioning or removing the sling from around the patient.

3. Lock the rear swivel casters.
4. Unhook the stand assist or transfer stand assist sling from all attachment points on the stand assist lift.
5. Instruct the patient to lift their feet off of the foot plate.

 Assist the patient if necessary.
6. Remove the stand assist or transfer stand assist from around the patient.
7. Unlock the rear casters and pull the stand assist lift away from the bed.

7 Troubleshooting

7.1 Troubleshooting Table

SYMPTOMS	FAULTS	SOLUTION
Noisy or dry sound from pivots.	Needs lubrication.	Refer to 8.6 Lubricating the Lift, page 40.
Electric actuator fails to lift or legs fail to open when button is pressed.	Hand-control or actuator connector loose.	Connect hand control or actuator connector. Ensure connectors are seated properly and fully connected.
	Battery low.	Charge batteries. Refer to 4.7 Charging the Battery, page 21.
	RED emergency stop button pressed IN.	Rotate RED emergency stop button CLOCKWISE until it pops out.
	Battery not connected properly to control box.	Reconnect the battery to the control box. Refer to 4.7 Charging the Battery, page 21.
	The connecting terminals are damaged.	Replace the battery pack. Refer to 4.7 Charging the Battery, page 21.
	Electric actuator in need of service or load is too high.	Refer to 8.7 Replacing the Electric Actuator, page 41 or 8.8 Replacing the Leg Actuators, page 42. Contact your Dealer or Invacare representative.
Unusual noise from lift arm actuator.	Actuator is worn or damaged or spindle is bent.	Refer to 8.7 Replacing the Electric Actuator, page 41 or 8.8 Replacing the Leg Actuators, page 42. Contact your Dealer or Invacare representative.
Lift arms will not lower in uppermost position.	Lift arms require a minimum weight load to lower from the uppermost position.	Pull down slightly on the lift arms.

SYMPTOMS	FAULTS	SOLUTION
Legs don't open and close properly.	Leg actuators may be worn or damaged or spindle is bent.	Refer to 8.8 Replacing the Leg Actuators, page 42.
Service light is flashing	12 months have passed since the lift was manufactured or initially assembled.	Initial assembly — Reset the service light. Refer to 3.5.1 Resetting the Service Light, page 16. After initial assembly — Service the lift. Contact your local Invacare dealer.

 If problems are not remedied by the suggested means, please contact your dealer or Invacare.

8 Maintenance

8.1 Maintenance Safety Inspection Checklist

WARNING!

- Maintenance **MUST** be performed **ONLY** by qualified personnel.
- A competent person must ensure that the strength and stability of the product remains adequate for all tasks if parts are replaced.
- **DO NOT** overtighten the mounting hardware, this will damage the mounting brackets.

Follow the maintenance procedures described in this manual to keep your patient lift in continuous service.

The Invacare Patient Lift is designed to provide a maximum of safe, efficient and satisfactory service with minimum care and maintenance.

All parts of the Invacare Lift are made of the best grades of steel, but metal to metal contact will wear after considerable use.

There is no adjustment or maintenance of the casters, other than cleaning, lubrication and checking axle and swivel bolts for tightness. Remove all debris, etc. from the wheel and swivel bearings. If any parts are worn, replace these parts **IMMEDIATELY**.

If you question the safety of any part of the lift, contact your Dealer **IMMEDIATELY**.

The UK Health and Safety Executive's Lifting Operations and Lifting Equipment Regulations 1998 require any equipment that is used in the workplace to lift a load be subject to safety inspection on a six monthly basis. This is to be carried out by a suitably trained individual. Please refer to the HSE website for guidance (www.hse.gov.uk).

Regular maintenance of patients lifts and accessories is necessary to assure proper operation. It is the responsibility of the person responsible for the patient lift to ensure that this maintenance is carried out.

After the first 12 months of operation, inspect all pivot points and fasteners for wear. If the metal is worn, the parts **MUST** be replaced. Perform this inspection every six months thereafter.

Date of Inspection:	Initials:
THE CASTER BASE <ul style="list-style-type: none"> <input type="checkbox"/> Inspect for missing hardware. <input type="checkbox"/> Base opens/closes with ease. <input type="checkbox"/> Inspect casters and axle bolts for tightness. <input type="checkbox"/> Inspect casters for smooth swivel and roll. <input type="checkbox"/> Inspect and clear wheels of debris. <input type="checkbox"/> Inspect pivot joints for wear. 	
SLINGS AND HARDWARE <ul style="list-style-type: none"> <input type="checkbox"/> Check all sling attachments each time it is used to ensure proper connection and patient safety. <input type="checkbox"/> Inspect sling material for wear. <input type="checkbox"/> Inspect straps for wear. <input type="checkbox"/> Inspect stitching. 	
THE ELECTRIC ACTUATOR ASSEMBLY <ul style="list-style-type: none"> <input type="checkbox"/> Check for leakage. <input type="checkbox"/> Inspect hardware on mast, boom and base. <input type="checkbox"/> Check for wear or deterioration. If damaged, return to factory. <input type="checkbox"/> Cycle to ensure smooth quiet operation of the electric actuator. 	

THE LIFT ARMS AND LINKAGE

- Check all hardware and attachment points.
- Inspect for bends or deflections.
- Inspect bolted joints of the lift arms for wear.
- Inspect to ensure that the lift arms are centered between the base legs.
- Inspect pivot joints for wear.
- Check sling hooks for wear or deflection.

THE MAST

- Mast must be securely assembled to lift arms.
- Inspect for bends or deflections.
- Inspect pivot joints for wear.

CLEANING

- Whenever necessary.

8.2 Service Interval

The lift should be serviced every 8000 cycles or annually, when the service light illuminates. Refer to 3.5 Checking the Service Light, page 16.

8.3 Cleaning the Sling and the Lift

WARNING!

- After EACH laundering (in accordance with laundering instructions on the sling), inspect the sling for wear, tears and loose stitching. Discard any patient sling that meets these criteria IMMEDIATELY.

Cleaning the Sling

The sling should be washed regularly in water temperature not exceeding 203°F (95°C) and a biological solution. Refer to the washing instructions in the sling user manual.

Cleaning and Disinfecting the Lift

To prevent cross-infection, the hoist must be cleaned and disinfected after each use.

A soft cloth, dampened with water and a small amount of mild detergent, is all that is needed to clean the patient lift. The lift can be cleaned with non-abrasive cleaners.

Never use acids, alkaline or solvents for cleaning the lift. Dry the lift carefully after cleaning.

Motors, control unit and mounting parts can be damaged if the lift is cleaned any other way than stated above.

The lift must be wiped with a moistened, firmly wrung cloth with ordinary household disinfectants. Only use disinfection detergents approved by the facility.

8.4 Reuse

This product is suitable for reuse. The maximum number of times it can be reused is dependent upon product condition.

To prevent the transmission of infection, the patient lift and slings must be cleaned after each use. Before reuse, refer to 8.3 Cleaning the Sling and the Lift, page 39.

8.5 Detecting Wear and Damage

It is important to inspect all stressed parts, such as slings, lifting arm and any pivot for slings for signs of cracking, fraying, deformation or deterioration. Replace any defective parts IMMEDIATELY and ensure that the lift is not used until repairs are made.


8.6 Lubricating the Lift

The Invacare lift is designed for minimum maintenance. However, a six month check and lubrication should ensure continued safety and reliability.

Keep lift and slings clean and in good working order. Any defect should be noted and reported to your dealer as soon as possible.

The casters MUST swivel and roll smoothly. A light grease (waterproof auto lubricant) may be applied to the ball bearing swivel of the casters once a year. Apply more frequently if the casters are exposed to extreme moist conditions.

8.7 Replacing the Electric Actuator


Ⓐ Electric Actuator	ⓓ Bearing
Ⓑ Actuator Lead Routing Hole	ⓔ Nylon Washer
Ⓒ Lift Arm Actuator Mounting Bracket	ⓖ Actuator Cap
Ⓓ Rubber Caps	ⓚ Bushing

Ⓔ Nut	Ⓛ Flat Washer
ⓕ Screw	Ⓜ Mast Assembly Actuator Mounting Bracket
Ⓖ Washer	

- Unplug the electric actuator from the control box on the mast assembly. Thread the actuator lead through the actuator lead routing hole in the mast assembly.
- Remove the hardware that secures the top of the electric actuator to the lift arm actuator mounting bracket.
 - Remove the rubber caps from the nut and screw that attach the electric actuator to the lift arm actuator mounting bracket.
 - Remove the screw, washer, bearing, nylon washer, actuator cap, bushing, nylon washer, bearing, washer, flat washer and nut.

 The actuator cap will become loose when the attaching hardware is removed. Inspect the actuator cap for damage. Replace the actuator cap if damage is found.
- Remove the hardware that secures the bottom of the electric actuator to the mast assembly actuator mounting bracket.
 - Remove the rubber caps from the nut and screw that attach the electric actuator to the mast assembly actuator mounting bracket.
 - Remove the screw, washer, bearing, nylon washer, bushing, nylon washer, bearing, washer, flat washer, and nut.
- To replace the electric actuator, reverse STEPS 1-3.

8.8 Replacing the Leg Actuators

Ⓐ Leg	Ⓓ Small Base Covers
Ⓑ Mast Handles	Ⓔ Pan Head Screws (eight screws attach the large base cover)
Ⓒ Screws (eight screws attach each small base cover)	Ⓕ Large Base Cover

If the legs of the patient lift do not open and close properly, it may be necessary to replace one or both of the leg actuators.

1. Test the operation of the legs:
 - Press the CLOSE button on the hand control to close the legs.
 - Press the OPEN button on the hand control to open the legs.
2. While holding on to the mast handles, carefully tip the patient lift back until it rests on the floor.

It is necessary to remove both small base covers before removing the large base cover.

3. Remove the sixteen screws that attach the small base covers to the base.
4. Remove the eight pan head screws that attach the large base cover to the base.

If the left or the right leg of the patient lift does not operate properly, perform the following steps on the side of the base necessary to replace the actuator. If both legs do not operate properly, perform this replacement procedure on both sides of the base.


5. Locate and disconnect the leg actuator lead wires.

The actuator lead wires will be visible after the large base cover is removed.

6. Remove the socket head screw and washer that attach the leg actuator to the base.
7. At the other end of the actuator, remove the socket head screw and nut that attach the leg actuator to the leg bracket.
8. If necessary, repeat STEPS 5-7 on the opposite side of the base to remove the other leg actuator.
9. To replace the actuator(s), reverse STEPS 5-7.

8.9 Adjusting the Knee Pad Height




A Knee Pad	B Quick Release Knob
------------	----------------------

i When adjusted properly, the patient's knee will make contact at the center of the knee pad.

1. While supporting the knee pad, disengage the quick release knob on the side of the knee pad mounting bracket.
2. Choose a height adjustment that will be comfortable for the patient and will provide the necessary support. Engage the quick release knob into position until the knee pad is secure in place.

9 Technical data

9.1 Patient Lift

Maximum height capacity (k)	161 cm
Lowest position (min. height) (l)	94 cm
Internal width at maximum reach (q)	86,5 cm
Total width (closed) external measure	67 cm
Base height (clearance)	12 cm
Base length	103,6 cm
Overall height	124,7 cm
Overall length	111,7 cm
Overall width	66 cm
Lifting area (height range) (m)	67 cm
Minimum internal width (r)	53,8 cm
Total width (open) internal measure	105 cm
Total width (open) centre to centre of castors (p)	109 cm
Turning radius	138 cm
Height to upper edge of legs	12 cm
Min. free height	4,5 cm
Maximum reach at 60 cm (a)	17,5 cm
Maximum reach from base (b)	41,4 cm
Reach from base with legs spread to 70 cm (c)	22,2 cm

Minimum Clearance of Base from the Floor (s)	3,7 cm
Maximum Clearance of Base from the Floor (t)	12 cm
Maximum External Length of the Base (u)	110,6 cm
Maximum Internal Length of the Base (v)	62,4 cm
Maximum External Width of Legs (y)	115,5 cm
Minimum Internal Width of Legs (z)	104,5 cm
Maximum Internal Width of Legs (w)	66,2 cm
Minimum Internal Width of Legs (x)	53,4 cm
Castor size (FRONT/REAR)	10 cm
Sling Options	Standing or Transfer
Sling Material	Polyester
Maximum weight limitation (SWL) (patient + sling + lift arms)	200 kg (31 stone)
Maximum patient weight	200 kg (31 stone)
Total weight (weight out of carton)	53 kg/117 lbs
Weight, Mast including battery	18,5 kg/40,8 lbs
Weight, base, leg section	20,5 kg/45,2 lbs
Weight, Knee Pad	7,5 kg/16,5 lbs

Weight, Foot Plate	6 kg/13 lbs
Operating forces of buttons	Max. 5N
Battery (voltage output)	24V DC mb max. 240VA (RCHBL)
Charger Input (voltage supply)	100-240V AC ~ 50/60 Hz
Charger Output/Charging Time	29.5V DC 2.9 Ah Max 6 hours
Audio/Visual Low Battery Alarm	Yes
Motor Safety Devices	Anti-Entrapment
*Approx. Lifts per Charge (working ability)	*100-200 Cycles per charge
Max. current input	Max. 400 mA
Operating temperature	5° to 40° C
Air humidity	20% to 90% at 30 not condensing
Sound pressure	< 50 dB
Atmospheric pressure	700 to 1060 hPa
Storage temperature	above 0° C
Storage air humidity	less than 60%
Storage atmospheric pressure	700 hPa to 1060 hPa
Degree of protection, control unit	IPX4
Degree of protection, hand control	IPX4
Degree of protection, motor	IPX4
Insulation class	Class II equipment, Type B applied part
Intermittens	10%, max, 2 minutes/18 minutes

Battery capacity	2,9 Ah
Manual emergency lowering	Yes
Electric emergency lowering/lifting	Yes/Yes
Expected lifetime	8 years
Lifting Speed	The Roze lift passes the “velocity of lifting and lowering” in the ENISO10535 (<0.15 m/s under maximum load and <0.25 m/s unloaded)

* Varies depending on load and stroke.

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Invacare Rehabilitation Equipment (Suzhou) Co. Ltd.
Invacare Corporation**USA**